Coagulation Calibrator Calibrateur de Coagulation

Code: HC00600 4 x 1 ml

Store at / Conservez entre 2 -8°C

Intended use

Coagulation calibrator is intended for calibration of haemostasis screening tests prothrombin time (PT; if to be expressed as %) and Fibrinogen. The values are given in the enclosed data sheet. For in vitro use and professional use only.

Summary

The PT test results can be reported in percentage, which means the proportional part of the normal PT activity (theoretical 100%), and is calculable from the calibration curve which plots the prothrombin time of different calibrator dilutions versus the inverse of the dilution.

The Fibrinogen concentration can be deduced from a log-log curve that plots the fibrinogen concentrations (g/l) of different calibrator dilutions versus the clotting time (sec).

It is possible to use the same calibration curve as long as the same lot of reagent is in use.

Reagents

It is a lyophilized, pooled human plasma with stabilizer.

Preparation

Rehydrate the vial with 1,0 ml of distilled water. Keep it at 18-25°C for 30 minutes. Swirl the vial gently before use and do not shake. Avoid contact of the fluid with the stopper.

Storage and stability

The intact vial is stable at 2 - 8°C up to the date of expiration as specified.

Stability after reconstitution: 4 hours at 22°C in the original vial. After further dilution, the diluted calibrator must be processed within 2 hours.

Procedure

Perform the calibration procedure as described in the insert for PT and Fibrinogen.

Precautions

- Do not use the reagent beyond the expiration date printed on the label.
- Avoid microbial contamination of the reagent or erroneous results may occur.
- Each donor unit used in the preparation of this reagent is tested and found to be negative for the following tests: antibodies to HIV, hepatitis C and hepatitis B surface antigen. However, the product must be handled with care, observing the precautions recommended for biohazardous material.

Bibliography

WHO Expert Committee on Standardisation: WHO Technical Report Series. No.889; 1999.

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