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Store at / Conservez entre 2 -8°C

### Intended use

The coagulation control (N+P) is intended for normal and pathological quality control for the following haemostasis screening tests PT, APTT and fibrinogen. It is not intended as a control for accuracy. For *in vitro* use and professional use only.

### Summary

The performance of daily quality control (QC) of haemostasis screening tests is highly recommended for security reasons. The industrial QC materials contain preservatives, stabilizers for lyophilization, which distort their original biological nature. For this reason they serve to check the repeatability and reproducibility of the measuring system. To check the accuracy of the reported result join and perform external quality assurance program in regular intervals.

### Reagents

It is a lyophilized, pooled human plasma with stabilizer. The pathological control is a modified form of the normal control, imitating pathological conditions of ill patients.

### Preparation

Rehydrate the vial with 1,0 ml of distilled water. Keep it at 18-25°C for 30 minutes. Swirl the vial gently before use and do not shake. Avoid contact of the fluid with the stopper.

### Storage and stability

The control, when stored at 2-8°C, will remain stable until the expiration date stated on the label. Once rehydrated, the reagent is stable for 4 hours at 22°C in the original vial.

### Procedure

The coagulation plasma must be tested in the same manner as any other plasma samples, according to the corresponding test instructions.

Compare the measured value with the declared values on the included sheet. The result has to be inside of the declared reference range.

### Precautions

- Each lot of these controls has been found to be negative for the following tests: antibodies to HIV, hepatitis C and hepatitis B surface antigen. However, the product must be handled with care, observing the precautions recommended for biohazardous material.
- Do not use the controls beyond expiration date printed on the label.
- Avoid microbiological contamination of the controls or erroneous results may occur.

### Expected values

The control range of the parameters may vary from lot to lot. The specific ranges are given on the included tables.

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