**GD-CAL****Multiparametric calibrator for clinical chemistry assays (human matrix)****REF**

GD8577 00

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GD-CAL is a human lyophilized serum. It is intended for the calibration of clinical chemistry methods with manual and automated procedures.

COMPOSITION**CALIBRATOR: 6 x 3 ml (lyophilized)**

Lyophilized human serum.
Biological additives and bacteriostatic agents.
The concentration and activities values of the components are indicated in the attached value sheet and are lot-specific.

PREPARATION

- Open the vial avoiding any loss of material
- Reconstitute with 3.0 ml of distilled water.
- Replace the stopper, gently swirl avoiding foam forming.
- Allow the calibrator to stand to room temperature for about 30 minutes before use.

Inaccurate reconstitution and/or storage and errors in assay technique can cause erroneous results.

STORAGE AND STABILITY

The calibrator is stable up to the expiration date stated on the label if stored tightly closed at 2-8 °C and contaminations are prevented during their use. Do not use the calibrator over the expiration date or if there is visible evidence of microbial growth. Store tightly capped when not use.

- After reconstitution is stable for:
 - 8 hours at 15-25 °C.
 - 2 days at 2-8 °C.
 - 1 month at -20 °C.
- Bilirubin (stored protected from light):
 - 4 hours at 15-25 °C.
 - 1 day at 2-8 °C.
 - 2 weeks at -20 °C.
- Acid Phosphatase:
 - 6 hours at +15-25 °C.
 - 2 days at +2-8 °C.
 - 2 weeks at -20 °C

Note:

1. Avoid freeze-thawing cycles.
2. It is suggested to subdivide reconstituted calibrator into aliquots and store at -20 °C.

PRECAUTIONS IN USE

All reagents contain inactive components such as preservatives (Sodium azide or others), surfactants etc. The total concentration of these components is lower than the limits reported by 67/548/EEC and 88/379/EEC directives about classification, packaging and labelling of dangerous substances. However, the reagents should be handled with caution, avoiding swallowing and contact with skin, eyes and mucous membranes.

The use of the laboratory reagents according to good laboratory practice is recommended.

The reagents from human donors have given negative results to anti-HIV 1/2, HBsAg and anti-HCV. It is recommended to handle with caution.

Waste Management

Please refer to local legal requirements.

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