# Polycheck<sup>®</sup>

**Scanner-based screening assay for the quantitative determination of allergen specific IgE**

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


## Materials Supplied, Storage and Stability

| Components  | Cat.-No.    | Content                         | Preparation                                     | Store at   | Shelf Life  |
|---|-------------|---------------------------------|---|--|---|
| <b>Polycheck® Pediatric Screen,</b><br>Chip coated with 5 calibrators and 20 allergens:<br>f13, f02, f01, f75, f35, f31, f03, f49, f14, f04, t03, g06, w06, d01, d02, e02, e01, e03, m03, m02     | MCPPC 1(Z)  | 20 (10) Chips                   | ready to use                                    | 2 - 8 °C<br><b>Protect from moisture!</b><br><br>Store together <b>with desiccant</b> in a carefully sealed <b>plastic bag</b> | until the expiration date                                 |
| <b>Polycheck® Inhalation Screen,</b><br>Chip coated with 5 calibrators and 20 allergens:<br>t03, t02, t04, t07, g06, g12, w06, w09, d01, d02, e02, e01, e03, e06, e84, e82, m03, m02, m01, m06    | MCIPC 1 (Z) | 20 (10) Chips                   |   |  |   |
| <b>Polycheck® Food Screen,</b><br>Chip coated with 5 calibrators and 20 allergens:<br>f17, f13, f16, f20, f02, f01, f75, f78, f35, f85, f31, f25, f03, f24, f95, f49, f14, f04, f10, f05          | MCFPC 1 (Z) | 20 (10) Chips                   |   |  |   |
| <b>Polycheck® Mediterranean Screen,</b><br>Chip coated with 5 calibrators and 20 allergens:<br>t03, t02, t04, t07, g06, g12, w01, w06, w09, w21, t09, d01, d02, e02, e01, e03, m03, m02, m01, m06 | MCMPC 1 (Z) | 20 (10) Chips                   |   |  |   |
| <b>Start Solution,</b><br>buffered protein solution   | MSPC        | 2 vials (1 vial)<br>5 mL, blue  | ready to use                                    | 2 - 8 °C   | 30 days after opening, respectively until expiration date |
| <b>anti-IgE Antibody,</b><br>monoclonal (murine) antibody labeled with ligand   | MAPC        | 2 vials (1 vial)<br>5 mL, green | ready to use                                    | 2 - 8 °C   | 30 days after opening, respectively until expiration date |
| <b>Enzyme-Labeled anti-Ligand,</b><br>Ligand conjugated to alkaline phosphatase   | MEPC        | 2 vials (1 vial)<br>5 mL, white | ready to use                                    | 2 - 8 °C   | 30 days after opening, respectively until expiration date |
| <b>Substrate Solution,</b><br>5'bromo-4'chloro-3' indolyl-phosphate-and 4' nitroblue tetrazolium, buffered  | MBPS3       | 2 vials (1 vial)<br>5 mL, black | ready to use                                    | 2 - 8 °C<br><b>protect from light!</b>   | 30 days after opening, respectively until expiration date |
| <b>Wash Buffer,</b><br>phosphate buffer, pH 7.4   | MWBPC       | 1 pouch                         | dissolve in 1 litre sterile demineralized water | 2 - 8 °C   | 30 days after opening, respectively until expiration date |

Material Safety Data Sheets are available on request (look as well [www.milenia-biotec.de](http://www.milenia-biotec.de)).

## Materials Required

- Personal computer (Windows 98, Windows ME, Windows 2000, Windows XP), printer
- Flat- bed- scanner (requested types on demand)
- Biocheck Imaging Software (BIS) for patient-oriented analysis, calculation and report
- shaker (orbital shaker: approx. 120 rpm)
- variable pipette for 250 µl up to 1000 µL; Multipette (1 mL) to dispense the wash solution
- sterile demineralized water
- sterile-one-litre bottle (for diluted wash solution)

## Specimen Collection and Preparation

Sera samples can be stored at 2-8 °C for up to 3 days. Storage at -20°C or lower is recommended for longer term storage. Avoid repeated freezing and thawing. Lipemic sera should be cleared by centrifugation.

## Warnings and Precautions

All reagents of this kit are strictly intended for *in-vitro* diagnostic use only. Do not replace or mix up reagents from one and the same or several test packages. Reagents' bottles shelf-life (30 days) starts with opening the bottles. Opened bottles should be used beforehand opening a new bottle! Only the powdered wash buffer can be exchanged between several kit packages. The application should be used by staff particularly informed and educated in the practice of *in-vitro* diagnostics. Compliance with the prescribed protocol for test implementation is always recommended. Sample material of patients (e.g. serum or plasma) normally used in laboratory determinations are always considered potentially infectious. Samples of such high risk patient groups should always be marked clearly and if necessary be handled in safety work benches (e.g. lamina flow bench).

## Method and Test Principle

The Polycheck<sup>®</sup> screening assay is an enzyme- immuno assay for the quantitative determination of allergen-specific IgE. There are 20 lines of relevant paediatric, food, or inhalant allergens coated together with 5 calibrators on a carrier membrane, which is located in the cavity of the Polycheck<sup>®</sup> chip. During incubation the patient's serum allergen-specific IgE binds to the corresponding allergens. Non-bound serum components are removed by washing. Ligand labeled anti-IgE antibodies bind to allergen-specific IgE. Unbound antibodies are removed by washing. Enzyme-labeled anti-ligand binds to the immune-complexes; surplus of enzyme conjugate is washed away. The substrate solution is added. Specifically bound enzymes convert the colourless substrate to a dark precipitate. The colour intensity of the (allergen-) lines is proportional to the respective allergen-specific IgE concentration in the patient's serum.

The Polycheck<sup>®</sup> chips can be evaluated using the Biocheck Imaging Software (BIS), a PC and a flat- bed scanner. The results are quantified in kU/L. This calculation is based on a number of calibrators, run on each chip.

## Test Performance

### Important notes:

- The required number of chips and test reagents should be at room temperature, before use.
  - Powdered wash buffer should be diluted with steril demineralized water at least 30 minutes prior to use. Use a steril bottle.
  - Do not allow the membrane of the test chips to dry during the assay is running.
  - All incubation steps are performed at room temperature (18-24°C) and with continual shaking.
1. Prepare a required number of Polycheck<sup>®</sup> allergen chips and mark them. Overlay allergen chips with **400 µL** Start Solution (blue) and incubate them for a few seconds. Tap the chips carefully on adsorbent paper.
  2. Add **400 µL** patient's serum and incubate for **45 minutes** on an orbital shaker.
  3. Decant sample, and wash **three times** with **1 ml** wash buffer.  
For that place the chips onto the work bench and pipet carefully 1 ml wash buffer into the cavity of the chip. Decant washing solution from the wells into a sink. Afterwards, add two times more 1 ml wash buffer in the same way into the cavities. After this, tap the chips carefully on adsorbent paper again.
  4. Then, pipet **400 µL** wash buffer into every chip cavity and incubate for **5 minutes** on the shaker. Decant and wash again every chip **3x with 1 ml** wash buffer. Tap the chips carefully on adsorbent paper.
  5. Pipet **400 µL** anti-IgE Antibody (green) and incubate for **45 minutes** on the shaker.
  6. Decant and wash every chip **3x with 1 ml** wash buffer. Tap carefully chips on adsorbent paper.
  7. Add **400 µL** Enzyme-labeled Anti-Ligand (white) and incubate for **20 minutes** on the shaker.
  8. Decant and wash every chip **3x with 1 ml** wash buffer. Tap carefully on adsorbent paper.
  9. Pipet **400 µL** substrate solution (black) and incubate **20 minutes in the dark** on the shaker.
  10. Decant and wash every chip **3x with 1 ml** wash buffer. Tap carefully the chips on adsorbent paper.
  11. Let the Polycheck<sup>®</sup>Chips become **air-dry**. Before evaluation with a scanner and the BIS software the membranes **must be absolutely dry!**

The chips can be evaluated up to 3 month, if they will be stored at room temperature and protected from light and humidity.

## Interpretation of Results

The evaluation of the Polycheck<sup>®</sup> allergen chips are operated with the assistance of a PC, the Biocheck Imaging Software (BIS), and a flat-bed scanner. The Software (BIS) analyzes, documents, and stores all patients' results in its database.

The Polycheck<sup>®</sup> Allergen chips are placed upside down onto the flat-bed-scanner's surface. Afterwards they are scanned and the data are processed by the software. - Please note orientation labelling on the surface.

The BIS-program verifies the technique of the chips, the validity of the calibrators, and ranges the intensity of greyscales and the position of the scanned band pattern to each allergen of the selected panel. BIS calculates allergen-specific IgE concentrations for each allergen on the membran. A patient's report is generated for every sample. There are all allergens listed which have been tested in this panel with the present concentration of allergen specific IgE. The results were quantitatively reported in International Units [kU IgE/L]. The RAST-class score to the respective IgE concentrations, and its clinical interpretation is shown in the following table:

| IgE [kU/L] | Class | Interpretation   |
|------------|-------|--|
| < 0.35     | 0     | No allergen specific antibody detectable   |
| 0.35 – 0.7 | 1     | Very low antibody titer; often no clinical symptom- despite sensitization        |
| 0.7 – 3.5  | 2     | Low antibody titer; sensitization, often clinical evidence in upper class ranges |
| 3.5 – 17.5 | 3     | Clear antibody titer; clinical evidence is mostly present                        |
| 17.5 – 50  | 4     | High antibody titer; always nearly with existing actuality clinical signs        |
| 50 – 100   | 5     | Very high antibody titer   |
| > 100      | 6     | Extreme high antibody titer  |

## Limitation of the Method

1. A definitive clinical diagnosis should not be based on the results of a single diagnostic method. It should only be made by a physician after all clinical and laboratory findings have been evaluated.
2. Allergenicity of food allergens may be modified by heating (cooking) and/or digestion. In this case the Polycheck<sup>®</sup> results may be inconsistent with the clinical findings.
3. Due to occurrence of common antigenic determinants in allergens of one natural substance family, the possibility of crossreactions cannot completely be excluded. With the Polycheck<sup>®</sup> system only the declared allergen panels can be covered, but not other (crossreacting) allergens!

## Assay Characteristics

|                        |  |
|------------------------|--|
| Sample material:       | serum  |
| Incubation time:       | 2 hours and 10 minutes at room temperature (18-24°C) on an orbital shaker            |
| Lower detection limit: | 0.15 kU/ L   |
| Specificity            | human IgE; no crossreaction with IgG, IgA, IgD, IgM                                  |
| Calibration:           | using single donor's and pool sera with defined allergen- specific IgE concentration |
| Method Comparison      | skin tests and alternative solid-phase immunoassays                                  |

## Short Instruction: Polycheck®

| Steps  | Solution                           | per Chip |
|--|------------------------------------|----------|
| Pipet  | Start Solution (blue)              | 400 µl   |
| Incubate a few seconds, decant and tap the chips on absorbant paper.   |                                    |          |
| Pipet  | Sample                             | 400 µl   |
| Incubate <b>45 min at room temperature (RT)</b> on an orbital shaker.  |                                    |          |
| Decant, wash <b>3x</b> with <b>1mL</b> wash buffer (per chip).<br>Add <b>400 µL</b> wash buffer and incubate <b>5 min</b> at <b>RT</b> on the shaker.<br>Decant, wash <b>3x</b> with <b>1mL</b> wash buffer, tap the chips on absorbant paper. | Wash Buffer                        | 400 µl   |
| Pipet  | Ligand-labeled anti-IgE (green)    | 400 µl   |
| Incubate <b>45 min at RT</b> on an orbital shaker.   |                                    |          |
| Decant, wash <b>3x</b> with <b>1mL</b> wash buffer (per chip), tap the chips on absorbant paper.   |                                    |          |
| Pipet  | Enzyme-labeled anti-Ligand (white) | 400 µl   |
| Incubate <b>20 min at RT</b> on a shaker.  |                                    |          |
| Decant, wash <b>3x</b> with <b>1mL</b> wash buffer (per chip), tap the chips on absorbant paper.   |                                    |          |
| Pipet  | Substrate Solution (black)         | 400 µl   |
| Incubate <b>20 min at RT</b> on a shaker <b>in the dark</b>  |                                    |          |
| Decant, wash <b>3x</b> with <b>1mL</b> wash buffer (per chip), tap the chips on absorbant paper.   |                                    |          |
| Air-dry the chips.   |                                    |          |
| Evaluate the <b>dry</b> chips by using Biocheck Imaging Software and a scanner.  |                                    |          |

For detailed description of the procedure see page 4, also.

## Available Panels:

| Pediatric (MKPPC)       |     | Inhalation (MKIPC)      |     | Food (MKFPC)         |     | Mediterranean (MKMPC)                         |     |
|-------------------------|-----|-------------------------|-----|----------------------|-----|---|-----|
| Peanut                  | f13 | Birch Pollen            | t03 | Hazelnut             | f17 | Birch Pollen / <i>Betula alba</i>             | t03 |
| Milk                    | f02 | Alder Pollen            | t02 | Peanut               | f13 | Alder Pollen / <i>Alnus spp.</i>              | t02 |
| Egg White               | f01 | Nutbush Pollen          | t04 | Walnut               | f16 | Nutbush Pollen / <i>Corylus avellana</i>      | t04 |
| Egg Yolk                | f75 | Oak Pollen              | t07 | Almond               | f20 | Oak Pollen / <i>Quercus robur</i>             | t07 |
| Potato                  | f35 | Timothy Grass           | g06 | Milk                 | f02 | Timothy Grass / <i>Phleum pratense</i>        | g06 |
| Carrot                  | f31 | Cultivated Rye          | g12 | Egg White            | f01 | Cultivated Rye / <i>Secale cereale</i>        | g12 |
| Codfish                 | f03 | Mugwort                 | w06 | Egg Yolk             | f75 | Common Ragweed / <i>Ambrosia elatior</i>      | w01 |
| Apple                   | f49 | English Plantain        | w09 | Casein               | f78 | Mugwort / <i>Artemisia vulgaris</i>           | w06 |
| Soybean                 | f14 | <i>D. pteronyssinus</i> | d01 | Potato               | f35 | English Plantain / <i>Plantago lanceolata</i> | w09 |
| Wheat Flour             | f04 | <i>D. farinae</i>       | d02 | Celery               | f85 | Pellitory / <i>Parietaria judaica</i>         | w21 |
| Birch Pollen            | t03 | Dog Epithelia           | e02 | Carrot               | f31 | Olive Pollen / <i>Olea europea</i>            | t09 |
| Timothy Grass           | g06 | Cat Epithelia           | e01 | Tomato               | f25 | <i>Dermatophagoides pteronyssinus</i>         | d01 |
| Mugwort                 | w06 | Horse Epithelia         | e03 | Codfish              | f03 | <i>Dermatophagoides farinae</i>               | d02 |
| <i>D. pteronyssinus</i> | d01 | Guinea pig Epithelia    | e06 | Shrimp               | f24 | Dog Epithelia                                 | e02 |
| <i>D. farinae</i>       | d02 | Hamster Epithelia       | e84 | Peach                | f95 | Cat Epithelia                                 | e01 |
| Dog Epithelia           | e02 | Rabbit Epithelia        | e82 | Apple                | f49 | Horse Epithelia                               | e03 |
| Cat Epithelia           | e01 | <i>A. fumigatus</i>     | m03 | Soybean              | f14 | <i>Aspergillus fumigatus</i>                  | m03 |
| Horse Epithelia         | e03 | <i>C. herbarum</i>      | m02 | Wheat Flour          | f04 | <i>Cladosporium herbarum</i>                  | m02 |
| <i>A. fumigatus</i>     | m03 | <i>P. notatum</i>       | m01 | Sesame               | f10 | <i>Penicillium notatum</i>                    | m01 |
| <i>C. herbarum</i>      | m02 | <i>A. tenius</i>        | m06 | Cultivated Rye Flour | f05 | <i>Alternaria tenuis</i>                      | m06 |

# Polycheck®

Scanner-based screening assay for the quantitative determination of allergen specific IgE  
English: Page 1-6

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