**BLOOD AGAR (Sheep Blood 7%)**

Medium for the cultivation and detection of haemolytic activity of streptococci and other fastidious microorganisms (ISO 10560: 1993).

**TYPICAL FORMULA (g/l)**

Proteose Peptone	15.0
Liver Digest	2.5
Sodium Chloride	5.0
Yeast Extract	5.0
Agar	12.0
Sheep Defibrinated Blood	70.0 ml
Final pH 7.2 ± 0.2	

**DESCRIPTION**

BLOOD AGAR (Sheep Blood 7%) is a medium for the cultivation and detection of haemolytic activity of streptococci and other fastidious microorganisms formulated according to ISO 10560: 1993.

**PRINCIPLE**

Proteose Peptone, liver digest and yeast extract provide essential carbon, nitrogen, vitamin and amino acid sources. Sodium chloride maintains the osmotic balance of the medium. Agar is the solidifying agent. The presence of blood provides further growth factors for fastidious microorganisms and allows to demonstrate the haemolytic reactions.

**TECHNIQUE**

Inoculate by streaking the specimen onto the surface of the medium in order to isolate pure cultures from specimens containing mixed flora. After streaking stab the agar several times to deposit beta-hemolytic streptococci beneath the agar surface. Incubate at 36±1°C for 18-48 hours, aerobically, anaerobically or under 5-10% CO<sub>2</sub> atmosphere, according to established laboratory procedures.

**INTERPRETATION of RESULTS**

Examine plates for growth and hemolytic reactions. Four types of hemolysis on blood agar media can be described:

1. α-hemolysis is the reduction of hemoglobin to methemoglobin in the medium surrounding the colony, causing a greenish discolorization of the medium.
2. β-hemolysis is the lysis of red blood cells, producing a clear zone surrounding the colony.
3. γ-hemolysis indicates no destruction of red blood cells and no change in the color of the medium.
4. δ-hemolysis indicates a partial lysis.

**STORAGE**

2-8°C away from light, until the expiry date on the label or until signs of deterioration or contamination are evident.

**WARNING and PRECAUTIONS**

The product is not classified as hazardous by current legislation and does not contain harmful substances in concentrations of ≥1%. The product is designed for *In vitro* diagnostic use and must be used only by properly trained operators.

**DISPOSAL of WASTE**

Disposal of waste must be carried out according to national and local regulations in force.

**REFERENCES**

1. ISO 10560: 1993. Milk and milk products – Detection of *Listeria monocytogenes*.
2. U.S. Food and Drug Administration. (1995). Bacteriological analytical manual, 8<sup>th</sup> ed., AOAC International, Gaithersburg, Md.

**Liofilchem s.r.l.**

Via Scozia-Zona industriale - 64026 Roseto degli Abruzzi Tel. +39.085.8930745 - Fax +39.085.8930330  
Web site: <http://www.liofilchem.net> E-mail: [liofilchem@liofilchem.net](mailto:liofilchem@liofilchem.net)

## PRODUCT SPECIFICATIONS

### NAME

BLOOD AGAR (Sheep Blood 7%).

### PRESENTATION

Ready plates (90 mm) containing 22+/-1 ml of medium.

### STORAGE

2-8°C

### PACKAGING

Code	Content	Packaging
10142	20 plates	<ul style="list-style-type: none"> <li>10 plates in thermically soldered film</li> <li>2 x 10 plates in cardboard boxes</li> </ul>
10142*	100 plates	<ul style="list-style-type: none"> <li>10 plates in thermically soldered film</li> <li>10 x 10 plates in cardboard box</li> </ul>

### pH OF THE MEDIUM

7.2 ± 0.2

### USE

BLOOD AGAR (Sheep Blood 7%) is a medium for the cultivation and detection of haemolytic activity of streptococci and other fastidious microorganisms formulated according to ISO 10560: 1993.

### TECHNIQUE

Refer to technical sheet of the product.

### APPEARANCE of the MEDIUM

Cherry red medium, opaque.

### SHELF LIFE

90 days.

### QUALITY CONTROL

- Control of general characteristics, label and print
- Sterility control
  - 7 days at 20 ± 2°C, in aerobiosis
  - 7 days at 35 ± 2°C, in aerobiosis
- Microbiological control
  - Inoculum for productivity: 10-100 UFC/ml
  - Inoculum for selectivity: 10<sup>4</sup>-10<sup>5</sup> UFC/ml
  - Inoculum for specificity: ≤ 10<sup>4</sup> UFC/ml
  - Incubation conditions: 18-24 h at 35 ± 2°C

Microorganism	ATCC	Growth	Characteristics
<i>Streptococcus pyogenes</i> *	ATCC 19615	good	β-hemolysis
<i>Streptococcus pneumoniae</i> *	ATCC 6303	good	α-hemolysis
<i>Neisseria meningitidis</i> *	ATCC 13090	good	γ-hemolysis
<i>Escherichia coli</i>	ATCC 25922	good	
<i>Staphylococcus aureus</i>	ATCC 25923	good	β-hemolysis

(\*in 10% CO<sub>2</sub> atmosphere)

### TABLE OF SYMBOLS

<b>IVD</b> <i>In vitro</i> Diagnostic Medical Device	<b>LOT</b> Batch code	 Manufacturer	 Contains sufficient for <n> tests
<b>REF</b> Catalogue number	 Temperature limitation	 Use by	 Caution, consult accompanying documents

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