COPRO SYSTEM

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COPRO SYSTEM

System for detection and presumptive identification of pathogenic intestinal microorganisms

DESCRIPTION

COPRO SYSTEM is an 18-well system containing desiccated biochemical substrates for the search for and presumptive identification of microorganisms from faecal samples.

The system is inoculated directly with faecal suspension and incubated at 36 ± 1 °C for 18-24 hours. The tests used in the search for and presumptive identification of the microorganisms present in the faecal sample are interpreted by assessing the colour change of the various wells and performing a microscope examination.

CONTENT OF THE PACKAGES

The packages contain:

- 20 COPRO SYSTEM (40 tests)
- 40 vials of *Inoculum Solution* (5.0 mL/vial)
- 1 sheet of instructions
- 1 TEST RESULTS FORM
- 20 Sealing films

ITEMS NECESSARY BUT NOT INCLUDED IN THE PACKAGES

- Oil of Vaseline for microbiological use code 80278
- OXIDASE TEST *STICK* code 88029
- Sterile containers for coproculture code 96708
- Cover slides
- Glass slides
- Microscope
- Physiological solution

CONFIGURATION

The configuration of the system is shown in Table 1. Each system allows examination of two samples.

Table 1

<i>Well</i>	BIOCHEMICAL REACTIONS FOR IDENTIFICATION
1-LDC	Decarboxylation of lysine
2-H₂S	Production of hydrogen sulphide
3-UR	Hydrolysis of urea
<i>Well</i>	PRESUMPTIVE IDENTIFICATION OF MICROORGANISMS
4-PRO	<i>Proteus spp./ Providencia spp.</i>
5-PSE	<i>Pseudomonas spp.</i>
6-YER	<i>Yersinia enterocolitica</i>
7-ESC	<i>Escherichia coli, E. coli O157</i>
8-KES	KES Group (<i>Klebsiella spp., Enterobacter spp., Serratia spp.</i>)
9-CAN	<i>Candida spp.</i>

PRINCIPLE OF THE METHOD

COPRO SYSTEM allows the search for and presumptive identification of microorganisms, pathogenic and otherwise, that may be present in faecal samples; for example: *Salmonella spp.*, *Citrobacter spp.*, *Proteus spp.*, *Pseudomonas spp.*, *Yersinia spp.*, *Escherichia coli*, *E.coli O157*, KES Group (*Klebsiella*, *Enterobacter*, *Serratia*), *Candida spp.*

- The presence of ***Salmonella spp.*** is shown by a change in colour from yellow to red of well **1-LDC**; a change in colour from yellow to black in well **2-H₂S** and by the yellow colour of well **3-UR**.
The identification of *Salmonella spp.* must be confirmed by isolation on selective media, biochemical tests and serological tests for typing of the group.
- The presence of ***Citrobacter spp.*** is shown by the yellow colour of well **1-LDC**; by the change of colour from yellow to black of well **2-H₂S** and by the yellow colour of well **3-UR**.
- The presence of ***Proteus spp./Providencia spp.*** is shown by the yellow colour of well **1-LDC**, by the change in colour from yellow to black of well **2-H₂S** and by the change in colour from yellow to red-fuchsia of well **3-UR**. Confirmation of *Proteus spp./Providencia spp.* is provided by the colour change from yellow to brown-black of well **4-PRO**.
- The presence of ***Pseudomonas spp.*** is shown by the change in colour from yellow to turbid green of well **5-PSE**. The identification of *Pseudomonas spp.* must be confirmed by the oxidase test (OXIDASE TEST STICK code 88029).
- The presence of ***Yersinia enterocolitica*** is shown by the change in colour from yellow to red-fuchsia of well **3-UR**, and by the change in colour from pink to turbid red of well **6-YER**. The identification of *Yersinia enterocolitica* must be confirmed by isolation on selective media, biochemical tests and serological typing tests.
- The presence of microorganisms of the group ***Escherichia coli, E. coli O157*** is shown by the change in colour from red to blue of the well **7-ESC**. The identification of *Escherichia coli* must be confirmed by a positive indole test, performed separately. The identification of *Escherichia coli O157* must be confirmed by isolation on selective media, biochemical tests and serological typing tests. The well also allows a search for *Shigella spp.*; in that case, a loop of broth culture is taken and sown on selective plates.
- The presence of microorganisms of the **KES GROUP** (*Klebsiella spp.*, *Enterobacter spp.*, *Serratia spp.*) is shown by the change in colour from violet to yellow of well **8-KES**. The identification of the microorganisms must be confirmed by isolation on selective media and by biochemical tests.
- The presence of ***Candida spp.*** is shown by the change in colour from green to yellow of the well **9-CAN**. The identification of *Candida spp.* must be confirmed by direct microscope examination of well 9, assessing the presence of mycelial hyphae and chlamydospores and by isolation on selective media.

COMPOSITION

Table 2

Well	Content
1-LDC	Culture medium with substrate for decarboxylation of lysine
2-H₂S	Culture medium with substrate for production of hydrogen sulphide
3-UR	Culture medium with substrate for hydrolysis of urea
4-PRO	Culture medium with substrate for isolation of <i>Proteus spp.</i> / <i>Providencia spp.</i>
5-PSE	Culture medium with substrate for isolation of <i>Pseudomonas spp.</i>
6-YER	Culture medium with substrate for isolation of <i>Yersinia enterocolitica</i>
7-ESC	Culture medium with substrate for isolation of <i>Escherichia coli, E. coli O157</i>
8-KES	Culture medium with substrate for isolation of the KES Group (<i>Klebsiella spp.</i> , <i>Enterobacter spp.</i> , <i>Serratia spp.</i>)
9-CAN	Culture medium with substrate for isolation of <i>Candida spp.</i>

Inoculum solution (g/L): Yeast extract **2.5 g**; Meat peptone **1.5 g**; Anhydrous glucose **1g**;
Distilled water **1000.0 mL**; pH **6.8 ± 0.2**

COLLECTION AND CONSERVATION OF THE SAMPLES

Obtain a sample with the specific sterile collection containers. In particular cases it may be useful to take the faecal material with a swab. The faecal samples to be subjected to the tests must be collected in accordance with the standard methods for samples to be submitted for coproculture.

TEST PROCEDURE

PREPARATION OF THE SAMPLE

Open one of the vials of Inoculum Solution contained in the kit.

1. FAECES

- Meticulously homogenise a sample of faeces of about 1 g in a test tube with 5 mL of physiological solution.
- Disperse 0.2 mL of homogenised sample in the vial of Inoculum Solution, shake and wait 5 minutes before sowing.

2. RECTAL SWAB

- Immerse the swab in a test tube containing 5 mL of physiological solution.
- Carefully squeeze the swab against the test tube wall so that the material to be examined is suspended homogeneously.
- Disperse 0.2 mL of sample in the vial of Inoculum Solution, shake and wait 5 minutes before sowing.

INOCULATION OF THE SYSTEM

1. Take a system from its wrapper and bring it to ambient temperature.
2. Write down the patient data and the date of the start of the examination.
3. Transfer 0.2 mL (4 drops) of suspension of the sample under examination into each well of the system.
4. Cover the first three wells: **1-LDC**, **2-H₂S**, **3-UR**, with three drops of oil of vaseline for microbiological use.
5. Cover the system with the lid provided and incubate at 36 ± 1 °C for 18-24 hours.

NB: The system is configured for performance of two tests. In case of performance of one test, cover the unused half system with the sealing film provided in the kit.

INTERPRETATION OF THE RESULTS

At the end of the incubation:

1. Take a drop of liquid from well **9-CAN**, deposit it on a glass slide, place a cover slide on top and examine under the microscope (40x) for the presence of chlamydo spores and hyphae.
2. Observe the changes of colour of the wells and interpret the results using table 3.
Note the results on the *TEST RESULTS FORM* (copy as many forms as necessary).

Table 3

Well	BIOCHEMICAL REACTIONS FOR IDENTIFICATION	Well colour		Microorganism			
		Positive reaction	Negative reaction	<i>Salmonella</i>	<i>Citrobacter</i>	<i>Proteus</i>	<i>Yersinia</i>
1-LDC	Decarboxylation of lysine	red +	yellow -	red +	yellow -	yellow -	yellow -
2-H₂S	Production of hydrogen sulphide	black +	yellow -	black +	black +	black +	yellow -
3-UR	Hydrolysis of urea	red-fuchsia +	yellow -	yellow -	yellow -	red+	red+

Well	PRESUMPTIVE IDENTIFICATION OF MICROORGANISMS	Well colour	
		Positive reaction	Negative reaction
4-PRO	<i>Proteus spp./ Providencia spp.</i>	brown-black	yellow
5-PSE	<i>Pseudomonas spp.</i>	turbid green	yellow-light blue
6-YER	<i>Yersinia enterocolitica</i>	turbid red	clear pink
7-ESC	<i>Escherichia coli, E.coli O157</i>	blue	grey-red
8-KES	<i>KES Group (Klebsiella spp., Enterobacter spp., Serratia spp.)</i>	yellow	violet
9-CAN	<i>Candida spp.</i>	yellow	green

+ : Positive reaction

- : Negative reaction

QUALITY CONTROL

Every batch of **COPRO SYSTEM** is subjected to quality control using the following reference microorganisms:

<i>Candida albicans</i>	ATCC 10231	<i>Proteus mirabilis</i>	ATCC 25933
<i>Citrobacter freundii</i>	ATCC 8090	<i>Pseudomonas aeruginosa</i>	ATCC 27853
<i>Escherichia coli</i>	ATCC 25922	<i>Salmonella typhimurium</i>	ATCC 14028
<i>Klebsiella pneumoniae</i>	ATCC 13883	<i>Yersinia enterocolitica</i>	ATCC 9610

FACTORS THAT MAY INVALIDATE THE RESULTS

- Poor standardisation of the inoculum.
- Material examined unsuitable.
- Use of expired systems.
- Non compliance with temperatures and times of incubation.

LIMITS AND WARNINGS

The presence of *Salmonella spp.*, *Escherichia coli*, *Yersinia spp.* and microorganisms of the KES Group (*Klebsiella spp.*, *Enterobacter spp.*, *Serratia spp.*) must be confirmed using appropriate culture media, biochemical tests (ENTEROSYSTEM 18R code 71618) and serological tests. For definitive identification of the microorganisms, it is necessary to make use of supplementary confirmation tests.

PERFORMANCE

The results obtained with **COPRO SYSTEM** agree with those obtained using traditional culture methods.

PRECAUTIONS

The product, **COPRO SYSTEM**, cannot be classified as hazardous under current legislation, nor does it contain harmful substances in concentrations $\geq 1\%$. It therefore does not require a Safety Data Sheet to be available. **COPRO SYSTEM** is a disposable device to be used only for diagnostic use *in vitro*. It must be used in the laboratory by properly trained operators, using approved aseptic and safety methods for handling pathogenic agents.

STORAGE

Store at 2-8 °C in the original packaging. Keep away from sources of heat and avoid excessive changes of temperature. In such conditions the product will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

DISPOSAL OF USED MATERIAL

After use, **COPRO SYSTEM** and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

PRESENTATION

Product	Code	Package
COPRO SYSTEM	71670	40 tests

TABLE OF SYMBOLS

 IVD <i>In Vitro</i> Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
 REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	 LOT Batch code

BIBLIOGRAFIA / BIBLIOGRAPHY

1. Murray, Baron, Pfaller, Tenorev and Tenover and Tenover *Manual of Clinical Microbiology 7th Edition (1999)* American Society of Microbiology
2. Ronald M. Atlas. *Handbook of Microbiological Media (1997)* Lawrence C. Parks

LIOFILCHEM S.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto D.A. (TE) - ITALY

Tel +390858930745 Fax +390858930330 Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net