# MULTODISCS

Anelli con otto antibiotici per test antibiogramma

## DESCRIZIONE

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# MULTODISCS

Rings with eight antibiotics for susceptibility tests

## DESCRIPTION

MULTODISCS are rings of absorbent paper with eight arms, each terminating in a 6 mm disc impregnated with an antibiotic. Each disc is marked with an identification code. A barrier of inert material at the end of each arm prevents diffusion of the antibiotic. This makes it possible to have easily measurable circular inhibition zones and makes MULTODISCS suitable for sensitivity tests using the NCCLS Bauer-Kirby method or other methods based on diffusion in agar media. MULTODISCS are offered in a large variety of configurations. Each configuration is available in the 100-test version.

## CONTENTS OF THE PACKAGES

Each pack contains: 100 rings packed in a jar with a dryer pack and an instruction leaflet.

## METHOD PRINCIPLE

The rings are applied to the surface of a culture medium inoculated with the culture broth prepared with pure colonies of the micro-organism under examination. After incubation the dishes are examined and the inhibition rings for each disc are measured and compared with the diameters of the standard inhibition ring. In this way the micro organism can be shown as sensitive, intermediate or resistant to the antimicrobial agents tested.

## COMPOSITION

The rings are made of high-quality paper in compliance with OMS and FDA specifications. Each disc is impregnated with a standardised amount of antibiotic.

## GATHERING AND KEEPING SAMPLES

The colonies that are to be subjected to the susceptibility testing are taken up by culture media that have been previously swabbed with the sample under examination. In the case of mixed colonies the bacterial strains must be purified before they are swabbed on the plates for the antibiogram.

## TEST PROCEDURE

For procedural details of the test, see the attached bibliography or the microbiology literature. A brief description of the NCCLS Bauer-Kirby method follows:

1. Take the jar from the refrigerator from the refrigerator and leave it for about an hour at ambient temperature in order to prevent condensation forming on the discs when it is opened as such condensation could affect long-term stability;
2. Swab 4 to 5 well isolated and morphologically similar colonies with a culture medium and suspend them in 5 ml of a suitable broth culture such as Tryptic Soy Broth;
3. Incubate the broth culture at  $36 \pm 1$  °C until 0.5 McFarland turbidity is reached;
4. Immerse a sterile buffer in the broth culture or in a diluted form thereof and squeeze it on the wall of the test tube to eliminate excess liquid. Drag it along the surface of the antibiogram medium contained on the plate so as to produce even growth;
5. Position the rings within 15 minutes from inoculation of the plates, pressing them with a sterile pliers on the surface of the agar and within 15 minutes of depositing of the discs leave the plate to incubate at  $36 \pm 1$  °C for 18-24 hours; incubate anaerobically at  $36 \pm 1$  °C for 48 hours, if anaerobic microorganisms are to be tested.

## EVALUATING THE RESULTS

At the end of the incubation period, measure the inhibition halos and compare them with the diameters of the standard NCCLS inhibition halos or with other standards used by the method adopted for the susceptibility testing.

## CLINICAL INTERPRETATION

The susceptibility testing carried out *in vitro* cannot exactly reproduce *in vivo* conditions. Nevertheless, it shows the effect of the concentration of the antibiotic, which varies in the culture medium in relation to the growth of the microbial population.

The final choice of antibiotic to administer to the patient is the responsibility of the clinician who possesses all the information on the patient.

## QUALITY CONTROL

Each batch of MULTODISCS is subjected to precise and thorough checks in compliance with NCCLS standards using the following bacterial strains:

<i>Bacteroides fragilis</i> ATCC 25285	<i>Haemophilus influenzae</i> ATCC 49247	<i>Pseudomonas aeruginosa</i> ATCC 27853
<i>Enterococcus faecalis</i> ATCC 29212	<i>Haemophilus influenzae</i> ATCC 49766	<i>Staphylococcus aureus</i> ATCC 25923
<i>Escherichia coli</i> ATCC 25922	<i>Klebsiella pneumoniae</i> ATCC 700603	<i>Streptococcus pneumoniae</i> ATCC 49619
<i>Escherichia coli</i> ATCC 35218	<i>Neisseria gonorrhoeae</i> ATCC 49226	

## LIMITS

Diffusion antibiograms use an *in vitro* technique and cannot therefore reproduce the extremely complex *in vivo* conditions. Nevertheless, it is a useful and important tool that helps the clinician choose the correct therapy. Many variable factors influence the final result of the diffusion antibiogram. The main ones are: the culture medium used, impregnation of the discs, inoculation of the medium, temperature, time and incubation atmosphere of the plates, pre-incubation and pre-diffusion conditions, depth of the medium, etc.

## PRECAUTIONS

The MULTODISCS cannot be classified as being hazardous according to current legislation but fall within the specific field of application where a safety data sheet must be supplied because they can cause phenomena of sensitization in sensitive subjects if they come into contact with the skin. MULTODISCS are disposable products. MULTODISCS are only for diagnostic *in vitro* use and are intended for professional use. They must be used in the laboratory by properly trained operators using approved aseptic and safety methods for pathogenic agents.

## STORAGE

MULTODISCS must be stored for long periods of time at -20 °C in their original packaging. For short periods, MULTODISCS can be stored at 2-8 °C, provided that they are still in their original packaging. Do not store them near sources of heat and do not expose them to excessive temperature variations. The packages must be placed in the refrigerator as soon as possible after use. In such conditions, MULTODISCS can be used until the expiry date shown on the label. Do not use after this date. Dispose of if they show signs of deterioration.

## ELIMINATING USED MATERIAL

After use, MULTODISCS and the material that comes into contact with the sample must be decontaminated and disposed of in accordance with current laboratory techniques for the decontamination and disposal of potentially infected material.

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