MYCOPLASMA SYSTEM *Plus*

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MYCOPLASMA SYSTEM *Plus*

System for the detection, semi-quantitative determination and susceptibility test of urogenital mycoplasmas

DESCRIPTION

MYCOPLASMA SYSTEM *Plus* is a 24-well system containing desiccated biochemical and antibiotic substrates for the detection, susceptibility test and semi-quantitative determination of urogenital mycoplasmas. The system is inoculated with the suspension of the clinical sample and incubated at 36 ± 1 °C for 18-24 hours. The tests used in the detection, count and susceptibility test are interpreted by evaluating the color of the various wells and performing a microscope examination.

KIT CONTENTS

The kit contains:

Ref. 72592	20 systems MYCOPLASMA SYSTEM <i>Plus</i>	20 Vials di <i>Physiological Solution</i> (7 mL/vial)
Ref. 79592	4 systems MYCOPLASMA SYSTEM <i>Plus</i>	4 Vials di <i>Physiological Solution</i> (7 mL/vial)
1 Instruction sheet		1 TEST RESULTS FORM

ITEMS NECESSARY BUT NOT INCLUDED IN THE KIT

Vaseline oil for microbiological use (VASELINE OIL 2 bottle of 50 mL) ref. 80278
Sundry microbiology laboratory materials

CONFIGURATION

The configuration of the system is shown in Table 1.

Table 1

Well	COUNT AND IDENTIFICATION OF MYCOPLASMAS / UREAPLASMAS
1-GR+	Growth of mycoplasmas ($10^2 < \text{titre} < 10^4$ CFU/mL)
2-GR++	Growth of mycoplasmas ($10^4 < \text{titre} < 10^5$ CFU/mL)
3-GR+++	Growth of mycoplasmas ($\text{titre} > 10^5$ CFU/mL)
4-ADC	Arginine Test
5-UR	Urea Test
Well	DETECTION OF T.VAGINALIS AND CANDIDA spp.
6-TR/YE	<i>Trichomonas vaginalis</i> / <i>Candida</i> spp.
Well	SUSCEPTIBILITY TESTING OF THE MYCOPLASMAS / UREAPLASMAS
7-TE	Tetracycline - 4 µg/mL
8-TE	Tetracycline - 8 µg/mL
9-PEF	Pefloxacin - 8 µg/mL
10-PEF	Pefloxacin - 16 µg/mL
11-OFX	Ofloxacin - 1 µg/mL
12-OFX	Ofloxacin - 4 µg/mL
13-DO	Doxycycline - 4 µg/mL
14-DO	Doxycycline - 8 µg/mL
15-E	Erythromycin - 8 µg/mL
16-E	Erythromycin - 16 µg/mL
17-CLA	Clarithromycin - 8 µg/mL
18-CLA	Clarithromycin - 16 µg/mL
19-MN	Minocycline - 4 µg/mL
20-MN	Minocycline - 8 µg/mL
21-CD	Clindamycin - 4 µg/mL
22-CD	Clindamycin - 8 µg/mL
23-AZM	Azithromycin - 4 µg/mL
24-AZM	Azithromycin - 8 µg/mL

PRINCIPLE OF THE METHOD

MYCOPLASMA SYSTEM Plus allows the detection, semi-quantitative count, presumptive identification and susceptibility test of *Mycoplasma hominis* and *Ureaplasma urealyticum* isolated from clinical samples and the detection and presumptive identification of the microorganisms most frequently isolated from vaginal and urethral swabs and seminal fluid, such as *Trichomonas vaginalis* and *Candida* spp.

- The semi-quantitative count of urogenital mycoplasmas (*Mycoplasma hominis* and *Ureaplasma urealyticum*) is provided by the yellow to red color change of wells **1-GR+**, **2-GR++**, **3-GR+++**.
- The presence of *Mycoplasma hominis* is shown by the yellow to red color change of the Arginine Test well **4-ADC**.
- The presence of *Ureaplasma urealyticum* is shown by the yellow to red color change of the Urea Test well **5-UR**.
- The presence of *Trichomonas vaginalis* and *Candida* spp. is determined by microscope (40x) examination of a drop of culture liquid taken from well **6-TR/YE** evaluating the presence of mobile ciliated protozoites for identification of *Trichomonas vaginalis* and/or the presence of chlamydospores and hyphae for identification of *Candida* spp.
- The susceptibility testing of the urogenital mycoplasmas is evaluated by observing the color change in wells **7-TE** to **24-AZM** as indicated in table 3.

COMPOSITION

Table 2

Well	Content
1-GR+	Culture medium for growth of mycoplasmas
2-GR++	Culture medium for growth of mycoplasmas
3-GR+++	Culture medium for growth of mycoplasmas
4-ADC	Culture medium with substrate for decarboxylation of arginine
5-UR	Culture medium con with substrate for hydrolysis of urea
6-TR/YE	Culture medium for isolation of <i>Trichomonas vaginalis</i> and <i>Candida</i> spp.
7-TE	Culture medium containing Tetracycline - 4 µg/mL
8-TE	Culture medium containing Tetracycline - 8 µg/mL
9-PEF	Culture medium containing Pefloxacin - 8 µg/mL
10-PEF	Culture medium containing Pefloxacin - 16 µg/mL
11-OFX	Culture medium containing Ofloxacin - 1 µg/mL
12-OFX	Culture medium containing Ofloxacin - 4 µg/mL
13-DO	Culture medium containing Doxycycline - 4 µg/mL
14-DO	Culture medium containing Doxycycline - 8 µg/mL
15-E	Culture medium containing Erythromicin - 8 µg/mL
16-E	Culture medium containing Erythromicin - 16 µg/mL
17-CLA	Culture medium containing Claritromycin - 8 µg/mL
18-CLA	Culture medium containing Claritromycin - 16 µg/mL
19-MN	Culture medium containing Minociclyne - 4 µg/mL
20-MN	Culture medium containing Minociclyne - 8 µg/mL
21-CD	Culture medium containing Clindamycin - 4 µg/mL
22-CD	Culture medium containing Clindamycin - 8 µg/mL
23-AZM	Culture medium containing Azithromycin - 4 µg/mL
24-AZM	Culture medium containing Azithromycin - 8 µg/mL

Physiological Solution (g/L): Sodium chloride **9 g**; Distilled water **1000 mL**; pH **6.8 ± 0.2**

COLLECTION AND CONSERVATION OF THE SAMPLES

Obtain the vaginal or urethral secretion with synthetic fibre swabs. Obtain the seminal fluid according to the technique envisaged for the microbiological culture examination. The samples must be sent for inoculation in the **MYCOPLASMA SYSTEM Plus** immediately after they are taken. Do not store in the refrigerator for even limited periods of time, since low temperatures may damage the vitality of some particularly sensitive microorganisms such as *Trichomonas vaginalis* prejudicing the final result.

TEST PROCEDURE**PREPARATION OF THE CLINICAL SAMPLE**

Open one of the vials of *Physiological Solution* contained in the kit.

- **VAGINAL SWAB - URETHRAL SWAB**

Immerse the swab (after obtaining the clinical material) in the vial of *Physiological Solution* and wait 5 minutes. Carefully squeeze the swab against the vial wall so that the clinical material is dispersed homogeneously in the physiological solution.

- **SEMINAL LIQUID**

Add 0.2 mL of sample to the vial of *Physiological Solution*, shake and wait 5 minutes before inoculating the system.

- **URINE**

Centrifuge 20 mL of urinary sample into two 10 mL tubes. Collect from each tube 0.1 mL of urinary sediment (0.2 mL) and add to the vial of *Physiological Solution* contained in the kit; shake and wait 5 minutes before inoculating the system.

- **CLINICAL MATERIAL FROM MYCOPLASMA TRANSPORT BROTH (ref. 20158)**

Add 1 mL of Mycoplasma Transport Broth, containing the sample under examination, to the vial of *Physiological Solution*, shake delicately and wait 5 minutes before inoculating.

INOCULATION OF THE SYSTEM

1. Take a system from its wrapper and allow to come to room temperature.
2. Write down the name of the patient, the date of the start of the examination, and the type of clinical material.
3. Transfer 0.2 mL (4 drops) of suspension of the clinical sample into each well of the system.
4. Cover all the wells except **6-TR/YE** well with one drop of Vaseline oil for microbiological use.
5. Cover the system with the lid provided and incubate at 36 ± 1 °C for 18-24 hours. If the presence of *Mycoplasma hominis* is suspected, incubate at 36 ± 1 °C for a further 24 hours.

INTERPRETATION OF THE RESULTS

At the end of the incubation:

- Watch for the color change of the wells **1-GR+** to **5-UR** and interpret the results using table 3 (COUNT AND IDENTIFICATION OF THE MYCOPLASMAS / UREAPLASMAS).
- Take a drop of liquid from well **6-TR/YE**, deposit it on a glass slide, place a cover slip on top and watch for the presence of *Trichomonas vaginalis* and *Candida* spp at the microscope (40x). Interpret as indicated in Table No. 3 (SEARCH FOR T. VAGINALIS AND CANDIDA spp).
- Watch for the color change of the wells from **6-TE** to **24-AZM** and interpret the results using Table No. 3 (SUSCEPTIBILITY TESTING OF THE MYCOPLASMAS / UREAPLASMAS). Refer to the well no. **1-GR+** as growth control for the antibiogram of mycoplasma.
- Note the results on the *TEST RESULTS FORM* (copy as many forms as necessary).

Table No.3

Well	COUNT AND IDENTIFICATION OF MICOPLASMAS / UREAPLASMAS	Well colour	
		Positive reaction	Negative reaction
1-GR+	Growth of mycoplasmas ($10^2 < \text{titre} < 10^4$ CFU/mL) ⁽¹⁾	red	yellow
2-GR++	Growth of mycoplasmas ($10^4 < \text{titre} < 10^5$ CFU/mL) ⁽²⁾	red	yellow
3-GR+++	Growth of mycoplasmas ($\text{titre} > 10^5$ CFU/mL) ⁽³⁾	red	yellow
4-ADC	Arginine Test: identifies <i>Mycoplasma hominis</i>	red	yellow
5-UR	Urea Test: identifies <i>Ureaplasma urealyticum</i>	red	yellow

Well	DETECTION OF <i>T. VAGINALIS</i> AND <i>CANDIDA</i> spp.	Microscopic observation (40x)
6-TR/YE	<i>Trichomonas vaginalis</i> / <i>Candida</i> spp.	<i>T. vaginalis</i> : mobile ciliated trophozoite <i>Candida</i> spp: chlamydo spores and hyphae

Well	SUSCEPTIBILITY TEST OF MICOPLASMAS / UREAPLASMAS	Well colour		
		S	I	R
7-TE	Tetracycline - 4 µg/mL	yellow	orange	red
8-TE	Tetracycline - 8 µg/mL	yellow	orange	red
9-PEF	Pefloxacin - 8 µg/mL	yellow	orange	red
10-PEF	Pefloxacin - 16 µg/mL	yellow	orange	red
11-OFX	Ofloxacin - 1 µg/mL	yellow	orange	red
12-OFX	Ofloxacin - 4 µg/mL	yellow	orange	red
13-DO	Doxycycline - 4 µg/mL	yellow	orange	red
14-DO	Doxycycline - 8 µg/mL	yellow	orange	red
15-E	Erythromycin - 8 µg/mL	yellow	orange	red
16-E	Erythromycin - 16 µg/mL	yellow	orange	red
17-CLA	Clarithromycin - 8 µg/mL	yellow	orange	red
18-CLA	Clarithromycin - 16 µg/mL	yellow	orange	red
19-MN	Minocycline - 4 µg/mL	yellow	orange	red
20-MN	Minocycline - 8 µg/mL	yellow	orange	red
21-CD	Clindamycin - 4 µg/mL	yellow	orange	red
22-CD	Clindamycin - 8 µg/mL	yellow	orange	red
23-AZM	Azithromycin - 4 µg/mL	yellow	orange	red
24-AZM	Azithromycin - 8 µg/mL	yellow	orange	red

S = Sensitive**I** = Intermediate sensitivity**R** = Resistant**LEGEND****(1)** = equivalent to 5-20 colonies grown on Mycoplasma Agar A7**(2)** = equivalent to 20-50 colonies grown on Mycoplasma Agar A7**(3)** = over 50 colonies grown on Mycoplasma Agar A7

QUALITY CONTROL

Every batch of **MYCOPLASMA SYSTEM Plus** is subjected to quality control using the following reference microorganisms:

<i>Mycoplasma hominis</i>	ATCC 23114
<i>Ureaplasma urealyticum</i>	ATCC 27618
<i>Candida albicans</i>	ATCC 10231

FACTORS THAT MAY INVALIDATE THE RESULTS

Poor standardisation of the inoculum; application of the method to microorganisms not in the Mycoplasma group; use of expired systems and/or reagents; incorrect application of the technique; non compliance with temperatures and times of incubation.

LIMITS AND WARNINGS

The quality of the collecting influences the test result. A negative test, therefore, is not always index of absence of infection. A positive test obtained with the system **MYCOPLASMA SYSTEM Plus** indicates an urogenital mycoplasmas infection, but it cannot be used, alone, to make a clinical diagnosis.

PERFORMANCE

The results obtained with **MYCOPLASMA SYSTEM Plus** agree with those obtained using traditional culture methods.

Application of the Wilcoxon non parametric statistical test shows that for each clinical material, the groups of standard values and those obtained with **MYCOPLASMA SYSTEM Plus** do not differ significantly $P \leq 0.05$ ⁽¹⁾.

PRECAUTIONS

The product **MYCOPLASMA SYSTEM Plus** is not classified as hazardous under current legislation, however see the safety data sheet for a correct use. **MYCOPLASMA SYSTEM Plus** is a disposable device to be used only for diagnostic use *in vitro*. It must be used in the laboratory by properly trained personnel, using approved aseptic and safety methods for handling pathogenic agents.

STORAGE

Store at 2-8 °C in the original packaging. Keep away from sources of heat and avoid excessive changes in temperature. In such conditions, the product will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

DISPOSAL OF USED MATERIAL

After use, **MYCOPLASMA SYSTEM Plus** and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

PRESENTATION

Product	Ref.	Packaging
MYCOPLASMA SYSTEM Plus	72592	20 tests
MYCOPLASMA SYSTEM Plus	79592	4 tests

TABLE OF SYMBOLS

IVD <i>In Vitro</i> Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	LOT Batch code

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LIOFILCHEM® S.r.l.