STAF SYSTEM 18R

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STAF SYSTEM 18R

System for the identification of *Staphylococcus* spp.

DESCRIPTION

STAF SYSTEM 18R is a 18 wells system containing desiccated biochemical substrata for the identification of *Staphylococcus* spp. The system is inoculated with the suspension of the microorganism to test and incubated at 36 ± 1 °C for 18-24 hours. The microorganism is identified by evaluating the color change in the different culture media in the wells and by a code number obtained from biochemical reactions interpretation.

CONTENT OF THE PACKAGES

The package contains:

20 Systems STAF SYSTEM 18R	20 Physiological Solution vials (7.0 mL/ vial)
1 Cartridge containing 20 xylose discs	20 Test Results Forms
1 Cartridge containing 20 ribose discs	1 Instruction sheet

NECESSARY ITEMS NOT CONTAINED IN THE PACKAGE

VASELINE OIL (ref. 80279)	NITRATE TEST (ref. 88009)
V.P. TEST – REAGENT (ref. 80280)	

CONFIGURATION

The system has the configuration indicated in table no. 1.

Table no.1

Well	BIOCHEMICAL REACTIONS FOR <i>Staphylococcus</i> spp. IDENTIFICATION
1- ADC	Arginine decarboxylase test
2- UR	Urea test
3- VP*	Indole test
4- NIT*	Nitrate test
5- ONPG	O-nitrophenyl-galattosidase test
6- NOV	Novobiocine test
7- MAL	Maltose test
8- TRE	Trehalose
9- MAN	Mannitol test
10- XYL	Xylose test
11- XYT	Xylitol test
12- CEL	Cellobiose test
13- SAC	Sucrose test
14- MNN	Mannose test
15- RIB	Ribose test
16- RAF	Raffinose test
17- LAC	Lactose test
18- FRU	Fructose test

: Add vaseline oil

* : Add the reagent indicated for the performance of the test

PRINCIPLE OF THE METHOD

STAF SYSTEM 18R allows the identification of *Staphylococcus* spp. of clinical significance. Identification is based on biochemical tests performed on specific culture media in the wells. The combination of positive and negative reactions allows to build up a code number which identifies bacteria by using the code number table.

COMPOSITION

Table no.2

Well	Content
1- ADC	Culture medium to detect arginine decarboxylation
2- UR	Culture medium to detect urea hydrolysis
3- VP	Culture medium to detect acetoin production
4- NIT	Culture medium to detect nitrates reduction
5- ONPG	Culture medium to detect O-nitrophenyl-galactopyranoside hydrolysis
6- NOV	Culture medium to detect resistance to novobiocine
7- MAL	Culture medium to detect maltose fermentation
8- TRE	Culture medium to detect trehalose fermentation
9- MAN	Culture medium to detect mannitol fermentation
10- XYL	Culture medium to detect xylose fermentation
11- XYT	Culture medium to detect xylitol fermentation
12- CEL	Culture medium to detect cellobiose fermentation
13- SAC	Culture medium to detect sucrose fermentation
14- MNN	Culture medium to detect mannose fermentation
15- RIB	Culture medium to detect ribose fermentation
16- RAF	Culture medium to detect raffinose fermentation
17- LAC	Culture medium to detect lactose fermentation
18- FRU	Culture medium to detect fructose fermentation

Physiological Solution (g/L): Sodium chloride 0,9 g; Distilled water 1000.0 mL; pH 7.0 ± 0.2

COLLECTION AND STORAGE OF THE SAMPLE

STAF SYSTEM 18R is used for the identification of bacteria belonging to the genus *Staphylococcus*, isolated on selective agar media such as Baird Parker Agar (ref. 10521), Vogel Johnson Agar (ref. 11185), Columbia C.N.A. Agar (ref. 11024), or non-selective media.

TEST PROCEDURE

PREPARATION OF THE BACTERIAL SUSPENSION

1. The microorganism to identify must have been recently isolated (18-24 hours); bacteria from more than 48 hours incubation can give not reliable results.
2. Open one vial of sterile physiological solution contained in the kit.
3. Take one or more morphologically similar, well-isolated colonies from the agar medium and suspend in the physiological solution so to obtain a turbidity equivalent to 1-2 McFarland (**Bacterial suspension**).
4. Thoroughly homogenize.

INOCULUM OF THE SYSTEM

1. Take a system from its pouch and leave it at room temperature.
2. Note date and origin of the microorganism.
3. Transfer 0.2 mL (4 drops) of **Bacterial suspension** into each well of the system and cover the wells **1-ADC** and **2-UR** with 2 drops of VASELINE OIL (ref. 80279).
4. Place a xylose disc in the well **10-XYL** and a ribose disc in the well **15-RIB**.
5. Cover the system with the lid provided and incubate at 36 ± 1 °C for 12-18-24 hours.

INTERPRETATION OF RESULTS

Subsequently to the incubation:

1. Introduce 2 drops of alpha-naphthol and 1 drop of NaOH 40% in the well **3-VP**. Wait for 15-20 minutes and watch for the appearance of a pink-red color.
2. Introduce 2 drops of sulphanic acid and 2 drops of alpha-naphthylamine in the well **4-NIT** (NITRATE TEST ref. 88009). Wait for a few seconds and watch for the appearance of a red-orange coloration.
3. Watch for color change of the wells and interpret the results using the table no. 3.
4. Note the result on the Test Results Form and build up the 6-digit code following the instructions provided in the paragraph **CODE NUMBER FORMING**.
5. Identify the microorganism using the code number table.

Table no. 3

Well	REACTIONS FOR BIOCHEMICAL IDENTIFICATION	Well color	
		Positive reaction	Negative reaction
1- ADC	Decarboxylation of arginine	Violet	Yellow
2- UR	Hydrolysis of urea	Fuchsia	Yellow
3- VP	Production of acetoin	Red- Pink	Yellow
4- NIT	Reduction of nitates	Red-Orange	Yellow
5- ONPG	Hydrolysis of ONPG	Yellow	Colorless
6- NOV	Resistance to novobiocine	Yellow- Orange	Red
7- MAL	Fermentation of maltose	Yellow- Orange	Red
8- TRE	Fermentation of trealose	Yellow- Orange	Red
9- MAN	Fermentation of mannitol	Yellow- Orange	Red
10- XYL	Fermentation of xylose	Yellow- Orange	Red
11- XYT	Fermentation of xylitol	Yellow- Orange	Red
12- CEL	Fermentation of cellobiose	Yellow- Orange	Red
13- SAC	Fermentation of sucrose	Yellow- Orange	Red
14- MNN	Fermentation of mannose	Yellow- Orange	Red
15- RIB	Fermentation of ribose	Yellow- Orange	Red
16- RAF	Fermentation of raffinose	Yellow- Orange	Red
17- LAC	Fermentation of lactose	Yellow- Orange	Red
18- FRU	Fermentation of fructose	Yellow- Orange	Red

CODE NUMBER FORMING

The biochemical tests are divided into 6 groups each containing 3 tests and each one is indicated with a positive value of 1, 2, 4.

- Value 1 : first test positive in each group (ADC, NIT, MAL, XYL, SAC, RAF).
- Value 2 : second test positive in each group (UR, ONPG, TRE, XYT, MNN, LAC).
- Value 4 : third test positive in each group (VP, NOV, MAN, CEL, RIB, FRU).
- Value 0 : every negative test

Adding the number of positive reactions in each group, it is obtainable a 6 digit code which, by the use of the code number table, allows the identification of the microorganism under examination as in the following example.

Example:

Test	GROUP 1			GROUP 2			GROUP 3			GROUP 4			GROUP 5			GROUP 6		
	ADC	UR	VP	NIT	ONPG	NOV	MAL	TRE	MAN	XYL	XYT	CEL	SAC	MNN	RIB	RAF	LAC	FRU
Positivity code	1	2	4	1	2	4	1	2	4	1	2	4	1	2	4	1	2	4
Results	-	+	+	+	-	-	+	+	+	-	-	-	+	+	+	-	-	+
Code	6			1			7			0			7			4		
CODE: 617074	IDENTIFICATION: <i>Staphylococcus aureus</i>																	

SIGNIFICANT BIOCHEMICAL REACTIONS FOR *Staphylococcus spp.* DIFFERENTIATION

	ADC	UR	VP	NIT	ONPG	NOV	MAL	TRE	MAN	XYL	XYT	CEL	SAC	MNN	RIB	RAF	LAC	FRU
<i>S. aureus</i>	V	+	+	+	V	-	+	+	+	-	-	-	+	+	+	-	V	+
<i>S. epidermidis</i>	V	+	+	V	V	-	+	-	-	-	-	-	+	V	V	-	V	+
<i>S. capitis</i>	V	-	V	V	-	-	-	-	+	-	-	-	V	+	-	-	-	+
<i>S. capae</i>	+	+	+	+	-	-	V	+	-	-	-	-	-	+	-	-	+	-
<i>S. warneri</i>	V	+	+	-	-	-	V	+	V	-	-	-	+	-	V	-	V	+
<i>S. haemolyticus</i>	+	-	V	V	-	-	+	+	V	-	-	-	+	-	V	-	V	V
<i>S. hominis</i>	V	+	V	V	-	-	+	V	-	-	-	-	V	-	-	-	V	+
<i>S. auricularis</i>	V	-	V	V	V	-	V	V	-	-	-	-	V	-	-	-	-	+
<i>S. saprophyticus</i>	-	V	+	-	V	V	+	+	V	-	V	-	+	-	-	-	V	+
<i>S. chonii var. human</i>	-	-	V	-	-	+	V	+	V	-	V	-	-	V	-	-	-	+
<i>S. chonii var. primate</i>	-	+	V	-	+	+	V	+	V	-	V	-	-	+	-	-	+	+
<i>S. xylosum</i>	-	+	V	V	+	+	+	+	V	+	V	-	+	+	V	-	V	+
<i>S. simulans</i>	+	+	-	+	+	-	-	V	+	-	-	-	+	V	-	-	+	+
<i>S. carnosus</i>	+	-	+	+	+	-	-	V	+	-	-	-	-	+	V	-	V	+
<i>S. intermedius</i>	V	+	-	+	+	-	V	+	V	-	-	-	+	+	+	-	V	+
<i>S. hyicus var. hyicus</i>	+	V	-	+	-	-	-	+	-	-	-	-	+	+	+	-	+	+
<i>S. hyicus var. chromogenes</i>	+	V	-	+	-	-	V	+	V	-	-	-	+	+	+	-	+	+
<i>S. sciuri</i>	-	-	-	+	-	+	V	+	+	-	-	+	+	V	+	-	V	+
<i>S. lentus</i>	-	-	-	+	-	+	V	+	+	V	-	+	+	V	+	+	V	V
<i>S. gallinarum</i>	-	+	-	+	-	+	+	+	+	+	V	+	+	+	+	+	V	+

+: Positive Reactions

-: Negative Reactions

V: Variable Reactions

QUALITY CONTROL

Each batch of **STAF SYSTEM 18R** is subjected to the quality control using the following reference microorganisms:

<i>Staphylococcus aureus</i>	ATCC 25923	<i>Staphylococcus xylosum</i>	ATCC 35033
<i>Staphylococcus sciuri</i>	ATCC 29062	<i>Staphylococcus epidermidis</i>	ATCC 12228

FACTORS THAT MAY INVALIDATE THE RESULTS

Poor standardization of the inoculum; unsuitable clinical material; use of expired systems and/ or reagents; non-compliance with incubation temperatures and times.

CAUTIONS

The product, **STAF SYSTEM 18R**, is not classified as hazardous under current legislation, however refer to the safety data sheet for a correct use. **STAF SYSTEM 18R** is a disposable device to be used only for *in vitro* diagnostic use. It must be used in the laboratory by properly trained personnel, using approved aseptic and safety methods for handling pathogenic agents.

STORAGE

Store **STAF SYSTEM 18R** at 2-8 °C in its original container. Keep away from sources of heat and avoid excessive changes of temperature. In such conditions **STAF SYSTEM 18R** maintains its validity until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

DISPOSAL OF USED MATERIAL

After use **STAF SYSTEM 18R** and the materials come into contact with the sample have to be decontaminated and disposed of in accordance with the techniques used in laboratory for decontamination and disposal of potentially infected material.

PRESENTATION

Product	Ref.	Package
STAF SYSTEM 18R	71630	20 tests

TABLE OF SYMBOLS

IVD <i>In vitro</i> Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contain sufficient for <n> tests	 Temperature limitation
REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	LOT Batch code

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LIOFILCHEM S.r.l.