

KJMR-IVA Blood Mixer and KJMR-IV Blood Mixer

KJMR-IVA Description:

KJMR-IVA is the upgrade model of KJMR-IV. Digitally controlling on the speed and time is applied to this new model on the basis of the old model function which is clearly displayed by LED. It would be more convenient for the operator and improve the accuracy of the experiment.

Specification:

Rotation Speed: 5-35 r/min

Timing range: 0-999min

Tube Qty: 1-24 pcs

Angle of Inclination: 38°

Voltage: AC 100-240V(50/60Hz)

Power Rating: 10W

Operation mode: Continuous/Timing



KJMR-IV Description:

The KJMR-IV Blood Mixer is specially designed for vacuum blood tube mixing. It is an ideal instrument of

blood corpuscle analyzer and blood flow analyzer for full sample mixing in laboratories with high mixing efficiency and consistency.

The instrument is made of silicone rubber which can well protect the surfaces of the vessels.

Specification:

Rolling Speed: 6-32 rpm

Tube Qty per : 1-24 pc

Angle of Inclination: 38°

Voltage: AC220V (50/60Hz)

Power Consumption: 12W

Working System: continue

Other Mixers are offered by Hasslab such as JJ-1 Series, KJMR series, QL Series and many other Magnetic mixers: 78-1 Series [Magnetic stir](#) and GL-3250 Series magnetic stirrer.

[Lab mixer](#) is the classic products of Hasslab, and magnetic stir is our hot sales, if you want to know more details about magnetic mixer stirrer, you can visit our Web or [contact us](#) as you like.