



MANUFACTURE STANDARD, MANUFACTURE RELIANCE

RT-150

Urine Analyzer



IDEAL CHOICE FOR

COMPACT SIZE

RELIABLE PERFORMANCE



Features

- Compact size, low weight.
- Cold light source, high luminosity with long life time.
- 5" Large LCD Display.
- Autosave data in case of power failure.
- Up to 1000 results memory.
- Build in thermal printer.
- Easy connection with PC by RS-232.
- Easy accessible design enables USER-DIY maintenance.



RT-150

Urine Analyzer



Test Report		
Sample No. 066		
11:13:24		2007-1-18
LEU	3+	500 cells/ul
BIL	3+	100 umol/L
BLD	3+	250 cells/ul
URO	Norm.	3.3 umol/L
KET	3+	8.0 mmol/L
GLU	3+	28 mmol/L
PRO	3+	3 g/L
PH	6	
NIT	2+	
SG	1.01	



Main Specifications

Test Items	Ascorbic Acid, Bilirubin, Blood, Glucose, Ketone, Leucocytes, Nitrite, PH, Protein, Specific gravity, Urobilinogen
Throughput	Normal mode, 60samples/hour Continuous mode, 120samples/hour
Display	5" Large LCD
Memory	1000 sample results
Output	Built-in thermal printer
Interface	RS-232 serial cable
Power supply	AC 110V - 220V ± 10%, 50-60Hz
Net Weight	2KG
Dimensions L × W × H (mm)	185 x 260 x 157

Strips for RT-150

11 test items, high sensitivity and specificity
Easy storage with special package against humidity

